



Top Pediatric Appliances

Available from models or scans

Fixed appliances available with doctor-fit bands, indirect bands, or SLS 3D Bands




Standard Expander

 Full palate scan




Slimline Expander

 Full palate scan




Nance

 Full palate scan




Lingual Arch

 Full palate scan




Band and Loop

 Quadrant scan




Habit Breaker

 Full palate scan




Sports Mouthguard

 Full palate scan





Theroux

 Full palate scan




3D Splints and Essix

 Horseshoe scan



KD Overlay

 Full palate scan

Top Scanning Advantages

● Improved patient experience

Say goodbye to the downsides of messy alignate! With digital impressions, you don't run the risk of warped impressions, long intraoral setting times, and you can start and stop your scan as needed to increase patient comfort.

● Indirect banding

Skip fitting bands chairside and save an appointment! Simply take and send your scan, and we'll fit bands on the teeth you designate.

● Combination orthodontics

Digital impressions can be especially useful for combination orthodontic treatments. We can remove brackets digitally, extract teeth as needed, or design whole appliances directly on the intraoral scan.

● Faster turnaround time

Eliminate inbound shipping time and get your cases into production, and back into your office, faster! Save even more time with our SLS 3D appliances.



Incoming Shipping

Up to 5 business days


Appliance Fabrication

Depending on appliance, usually 10 business days in lab

Return Shipping

1-2 days


Bands Options



DOCTOR FIT

Fit your bands chairside and send them in to the lab. We'll seat on your model.

IMPRESSION: With bands
SCAN: Without bands



INDIRECT

Skip the fitting, and let us fit and place bands for you on the teeth you specify.

IMPRESSION: Without bands
SCAN: Without bands




3D / SLS

We'll design custom CAD bands directly on your intraoral scan and print in stainless steel!

IMPRESSION: *Not available.*
SCAN: Without bands

Sending Prescriptions



NEOLab Junior is home to RX Wizard, a proprietary prescription management system. Here you can submit new prescriptions, save your favorite appliance designs to use again and again, and batch your cases to save on shipping.

● Templates

Save your favorite appliance designs! Just add a patient name, due date, and submit.

● Batch Shipping

See when return shipments are scheduled and add as you like to them. Shipping is by box, not case, so fill your box as full as you can!

● Case Transparency

Constant status updates on your case mean you can see when it's arrived, when it's in production, and when it's on the way back to your office - with integrated UPS tracking!